



503

3 Section Plinth

Very popular and versatile three-section couch, firm and stable providing the ideal base for general treatment techniques through to heavy manipulation.

Features: unique lifetime guarantee, easy operation, hoist accessible, low running costs, high density fire retardant foam, hardwearing anti-microbial vinyl, able to lift 225kg from its lowest height, full choice of accessories and colours.

STANDARD SPECIFICATION

- Either electric or hydraulic height adjustment available.
- Able to lift up to 225kg (35 stone) from its lowest height.
- Backrest & footrest adjustment by gas strut mechanism.
- Backrest & footrest angle adjustable from -30° to + 85° vertical.
- Hoist clearance of 13cm (5").
- Adjustable foot for uneven floors.
- Fully washable/ fire retardant Ambla vinyl available in a range of colours.
- Couch goes horizontal for patient recovery in emergency situations.
- Robust retractable castors for excellent mobility and safety.
- CE marked in conformity with MDD 93/42/EEC Annex IX Class 1.
- Electrical standard IEC601 part 1.


STANDARD OPTIONS

- Breathing hole and plug.
- Motorised backrest (electric versions only).
- Battery backup (electric versions only).
- Cot Sides.
- Paper roll holder.
- Extra width upholstery up to 74cm (29").
- Double thickness padded upholstery.
- Foot switch.
- Foot ring switch.
- Other options available.

503 DIMENSIONS

Length: 185 cm (73")
 Width: 63 cm (25")
 Height: 46 cm (18") approx
 to 95 cm (37") approx

503 Electric
 503 Hydraulic



Model shown with optional breathing hole.

PACKED WEIGHTS AND DIMENSIONS FOR SHIPPING

Packed Weight: 82 Kgs
 Shipping Dimensions: 190 x 66 x 64 cm

DELIVERY PROMISE

Our normal delivery service is 7 – 10 days from receipt of order; however if you have a specific requirement we will always endeavour to meet it.

FOR YOUR PEACE OF MIND

Our unique guarantee is for a LIFETIME on the mechanism and UNCONDITIONAL for 3 years on all components.. This even applies to accidental damage. We will not quibble. The only exclusion is damage to upholstery.

All of our products are CE marked. We are an ISO9001 audited company. Should you require any help concerning the most suitable product for your needs, we would be delighted to discuss it with you.



Plinth 2000 Ltd.

Wetheringsett Manor
 Wetheringsett
 Stowmarket
 Suffolk IP14 5PP
 United Kingdom

t: +44 (0)1449 767887
 f: +44 (0)1449 766122
 e: sales@plinth2000.com

www.plinth2000.com

Plinth | 2000



VINYL

Upholstery

We use only the best quality vinyl proven over many years and **ALL** feature an anti microbial additive to help prevent cross infection.

All coverings are vinyl and washable, being fire retardant to BS 7176, BS 5852 part 2 (crib 5) and display excellent levels of wear. ALL include patented anti microbial technology proved effective in inhibiting the growth of E-Coli and Staphylococcus Aureus (MRSA), and all are suitable for cleaning and disinfecting with NaDCC based products.

Note: Unless specified BEIGE will be supplied.



Beige




Grey




Black



Fawn



Chocolate




Brow




Damson




Wine



Ruby



Terracotta



Night




Ocean




Sky



Slate




Green



Sage



Lime



White

Plinth 2000 Ltd.

Wetheringsett Manor
Wetheringsett
Stowmarket
Suffolk IP14 5PP
United Kingdom

t: +44 (0)1449 767887
f: +44 (0)1449 766122
e: sales@plinth2000.com

www.plinth2000.com

Plinth | 2000